

For Research Use Only

IQGAP2 Polyclonal antibody, PBS Only

Catalog Number:30692-1-PBS



Basic Information

Catalog Number:

30692-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33518

GenBank Accession Number:

NM_006633

GeneID (NCBI):

10788

UNIPROT ID:

Q13576

Full Name:

IQ motif containing GTPase activating protein 2

Calculated MW:

181 kDa

Observed MW:

181 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

IQGAP2, a member of the IQGAP family, functions as a tumor suppressor in most of the cancers. IQGAP2 contains three IQ domains, one calponin homology domain, one Ras-GAP domain and one WW domain. IQGAP2 interacts with components of the cytoskeleton, with cell adhesion molecules, and with several signaling molecules to regulate cell morphology and motility. It also acts as a tumor suppressor and has been found to play a role in regulating innate antiviral responses.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30692-1-AP (IQGAP2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30692-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 30692-1-AP (IQGAP2 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 30692-1-PBS in a different storage buffer formulation.

Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 30692-1-AP (IQGAP2 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 30692-1-PBS in a different storage buffer formulation.