

For Research Use Only

TCF25 Polyclonal antibody

Catalog Number: 30852-1-AP



Basic Information

Catalog Number:

30852-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29375

GenBank Accession Number:

BC000959

GeneID (NCBI):

22980

UNIPROT ID:

Q9BQ70

Full Name:

transcription factor 25 (basic helix-loop-helix)

Calculated MW:

77 kDa

Observed MW:

70-80 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:2000-1:12000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : rat heart tissue, mouse heart tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30852-1-AP (TCF25 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.