

For Research Use Only

TNIK Polyclonal antibody, PBS Only

Catalog Number: 30906-1-PBS



Basic Information

Catalog Number: 30906-1-PBS	GenBank Accession Number: BC150256	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 23043	
Source: Rabbit	UNIPROT ID: Q9UKE5	
Isotype: IgG	Full Name: TRAF2 and NCK interacting kinase	
Immunogen Catalog Number: AG34132	Observed MW: 155-180 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

TNIK (Traf2- and Nck-interacting kinase) is one of the germinal center kinase family members involved in cytoskeleton organization and neuronal dendrite extension, which had been implicated in postsynaptic signaling as well as in regulation of cell proliferation. TNiK is expressed in the nervous system and plays key roles at the synapse and nucleus. In non-neuronal cells, TNiK has been described as a critical component of transcriptional regulatory mechanisms, mediated by Wnt signaling (PMID: 24566388). TNiK is also an essential regulatory component of the T-cell factor-4 and β -catenin transcriptional complex and is required for the tumor-initiating function of colorectal cancer stem cells (PMID: 10521462). Western blot analysis detected TNiK at an apparent molecular mass of 150-180 kDa as cited in literature (PMID: 24566388).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30906-1-AP (TNIK antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30906-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed K-562 cells using TNIK antibody (30906-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 30906-1-PBS in a different storage buffer formulation.