

For Research Use Only

AHCTF1 Polyclonal antibody

Catalog Number: 30923-1-AP



Basic Information

Catalog Number:

30923-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32788

GenBank Accession Number:

NM_001323342.1

GeneID (NCBI):

25909

UNIPROT ID:

Q8WYP5

Full Name:

AT hook containing transcription factor 1

Calculated MW:

252 kDa

Observed MW:

270 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, A549 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30923-1-AP (AHCTF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.