

For Research Use Only

AHCTF1 Polyclonal antibody

Catalog Number: 30923-1-AP **1 Publications**



Basic Information

Catalog Number: 30923-1-AP	GenBank Accession Number: NM_001323342.1	Purification Method: Antigen affinity purification
Source: Rabbit	GeneID (NCBI): 25909	Recommended Dilutions: WB: 1:5000-1:50000 IF/ICC: 1:400-1:1600
Isotype: IgG	UNIPROT ID: Q8WYP5	
Immunogen Catalog Number: AG32788	Full Name: AT hook containing transcription factor 1	
	Calculated MW: 252 kDa	
	Observed MW: 270 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls: WB : HeLa cells, A549 cells
Cited Applications: WB	IF/ICC : A431 cells,
Species Specificity: human	
Cited Species: human	

Background Information

AHCTF1, also known as ELYS, TMBS62, is a large multifunctional protein. It plays essential roles in nuclear pore assembly and mitosis. AHCTF1 was expressed at higher levels in embryonic liver, spleen and thymus compared to the corresponding adult tissues (PMID: 11952839) and has 3 isoforms with MW of 252-256 kDa.

Notable Publications

Author	Pubmed ID	Journal	Application
Huifang Liu	40707755	Dig Dis Sci	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30923-1-AP (AHCTF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using AHCTF1 antibody (30923-1-AP) at dilution of 1:800 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002), CL594-phalloidin (red).