

For Research Use Only

AMAC1L2 Polyclonal antibody, PBS Only

Catalog Number:30979-1-PBS



Basic Information

Catalog Number:

30979-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34353

GenBank Accession Number:

BC131696

GeneID (NCBI):

83650

UNIPROT ID:

Q96KT7

Full Name:

acyl-malonyl condensing enzyme 1-like 2

Calculated MW:

388 aa, 35 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

AMAC1L2, also known as SLC35G5, belongs to the transporter protein family 35 member G5. This protein has an N-terminal α -helix, a central β -side and a C-terminal structural domain at the cell membrane. The role of SLC35G5 in the tumor microenvironment as a drug target and also as a biomarker for diagnosis and treatment monitoring.

Storage

Storage:

Store at -80°C .

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30979-1-AP (AMAC1L2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30979-1-PBS in a different storage buffer formulation.