

For Research Use Only

FPR2 Polyclonal antibody, PBS Only

Catalog Number: 30989-1-PBS



Basic Information

Catalog Number: 30989-1-PBS	GenBank Accession Number: BC029125	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 2358	
Source: Rabbit	UNIPROT ID: P25090	
Isotype: IgG	Full Name: formyl peptide receptor 2	
Immunogen Catalog Number: AG34304	Calculated MW: 351 aa, 39 kDa	
	Observed MW: 42 kDa	

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

The formyl peptide receptors (FPRs) are a group of G protein-coupled chemoattractant receptors that play important roles in host defense and inflammation (PMID: 28335409). In humans, three different isoforms are expressed (FPR1, FPR2, and FPR3). FPR2, also known as FPRL1 or ALX, is expressed by many cell types including monocytes, neutrophils, macrophages, astrocytoma cells, and epithelial cells. FPR2 conveys the biological functions of a variety of ligands, including N-formyl-methionyl peptides, the pro-resolution mediators annexin A1 (AnxA1) and lipoxin A4, as well as the activating and proinflammatory protein serum amyloid A (PMID: 22610094).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 30989-1-AP (FPR2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30989-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 30989-1-AP (FPR2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 30989-1-PBS in a different storage buffer formulation.