

For Research Use Only

PCYT1A Polyclonal antibody, PBS Only

Catalog Number: 30991-1-PBS



Basic Information

Catalog Number: 30991-1-PBS	GenBank Accession Number: NM_005017	Purification Method: Antigen affinity Purification
Source: Rabbit	GeneID (NCBI): 5130	
Isotype: IgG	UNIPROT ID: P49585	
Immunogen Catalog Number: AG34268	Full Name: phosphate cytidyltransferase 1, choline, alpha	
	Calculated MW: 42 kDa	
	Observed MW: 40-41 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

PCYT1A (phosphate cytidyltransferase 1A, choline), also known as CTPCT. It is expected to be located in the nucleus and cytoplasm, which is enriched in brain, placenta, liver, fetal and adult lung. This gene belongs to the cytidyltransferase family and is involved in the regulation of phosphatidylcholine biosynthesis. Mutations in this gene are associated with spondylometaphyseal dysplasia with cone-rod dystrophy. The molecular weight of PCYT1A is 41 kDa.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 30991-1-AP (PCYT1A antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30991-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 30991-1-AP (PCYT1A antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30991-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using PCYT1A antibody (30991-1-AP) at dilution of 1:400 and Coralite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 30991-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed U-251 cells using PCYT1A antibody (30991-1-AP) at dilution of 1:400 and Coralite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 30991-1-PBS in a different storage buffer formulation.