

For Research Use Only

# MFRP Polyclonal antibody, PBS Only

Catalog Number: 31017-1-PBS



## Basic Information

Catalog Number:

31017-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34302

GenBank Accession Number:

BC160191

GeneID (NCBI):

83552

UNIPROT ID:

Q9BY79

Full Name:

membrane frizzled-related protein

Calculated MW:

62 kDa

Observed MW:

62 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31017-1-AP (MFRP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31017-1-PBS in a different storage buffer formulation.