

For Research Use Only

Collagen Type XV Polyclonal antibody, PBS Only

Catalog Number:31056-1-PBS



Basic Information

Catalog Number:

31056-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34644

GenBank Accession Number:

BC157094

GeneID (NCBI):

1306

UNIPROT ID:

P39059

Full Name:

collagen, type XV, alpha 1

Observed MW:

250 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31056-1-AP (Collagen Type XV antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31056-1-PBS in a different storage buffer formulation.