

For Research Use Only

GIGYF1 Polyclonal antibody

Catalog Number: 31073-1-AP



Basic Information

Catalog Number:

31073-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34632

GenBank Accession Number:

BC129991

GeneID (NCBI):

64599

UNIPROT ID:

O75420

Full Name:

GRB10 interacting GYF protein 1

Observed MW:

140 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB: HEK-293T cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293T cells were subjected to SDS PAGE followed by western blot with 31073-1-AP (GIGYF1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.