

For Research Use Only

GIGYF1 Polyclonal antibody, PBS Only

Catalog Number: 31073-1-PBS



Basic Information

Catalog Number: 31073-1-PBS	GenBank Accession Number: BC129991	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 64599	
Source: Rabbit	UNIPROT ID: O75420	
Isotype: IgG	Full Name: GRB10 interacting GYF protein 1	
Immunogen Catalog Number: AG34632	Observed MW: 140 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293T cells were subjected to SDS PAGE followed by western blot with 31073-1-AP (GIGYF1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31073-1-PBS in a different storage buffer formulation.