

For Research Use Only

ADAMTS14 Polyclonal antibody, PBS Only

Catalog Number: 31075-1-PBS



Basic Information

Catalog Number:

31075-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34638

GenBank Accession Number:

NM_080722

GeneID (NCBI):

140766

UNIPROT ID:

Q8WXS8

Full Name:

ADAM metalloproteinase with thrombospondin type 1 motif, 14

Calculated MW:

134 kDa

Observed MW:

134 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

The ADAMTS (a disintegrin-like and metalloproteinase with thrombospondin motifs) is a family of extracellular metalloproteinases mediating diverse functions including matrix degradation, blood coagulation, and angiogenesis. The ADAMTS family constitutes a group of proteins composed of 19 enzymes and 7 ADAMTS-like proteins. ADAMTS14 possesses high similarity with ADAMTS2 and ADAMTS3 in sequence and domain composition (PMID: 11741898). The expression of ADAMTS14 is correlated with the occurrence and progression of colorectal cancer and oral squamous cell carcinoma (PMID: 34856162; PMID: 32102222).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31075-1-AP (ADAMTS14 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31075-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 31075-1-AP (ADAMTS14 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31075-1-PBS in a different storage buffer formulation.