

For Research Use Only

ANK2 Polyclonal antibody, PBS Only

Catalog Number: 31080-1-PBS



Basic Information

Catalog Number:

31080-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34642

GenBank Accession Number:

NM_001148

GeneID (NCBI):

287

UNIPROT ID:

Q01484

Full Name:

ankyrin 2, neuronal

Calculated MW:

434 kDa

Observed MW:

250-260 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IHC, IF-P, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

The Ankyrins, belonging to the adapter protein family, consist of Ankyrin-R (ANK1), Ankyrin-B (ANK2), and Ankyrin-G (ANK3). Initially, ANK2 was considered as a key cytoskeletal-associated protein in the brain. Besides, it was discovered as a pivotal adaptor and scaffolding protein in the heart. Remarkably, recent studies have demonstrated that Ankyrins contribute to the modulation of various biological activities in cancers. (PMID: 36539782)

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 31080-1-AP (ANK2 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31080-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse cerebellum tissue slide using 31080-1-AP (ANK2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31080-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse heart tissue slide using 31080-1-AP (ANK2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31080-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse heart tissue using ANK2 antibody (31080-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31080-1-PBS in a different storage buffer formulation.