

For Research Use Only

PYROXD1 Polyclonal antibody, PBS Only

Catalog Number: 31093-1-PBS

Featured Product



Basic Information

Catalog Number:

31093-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34768

GenBank Accession Number:

BC021662

GeneID (NCBI):

79912

UNIPROT ID:

Q8WU10

Full Name:

pyridine nucleotide-disulphide oxidoreductase domain 1

Observed MW:

55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

pyridine nucleotide-disulfide oxidoreductase domain 1 (PYROXD1) is a putative oxidoreductase that co-evolved with the tRNA-LC and is essential for its activity in cells and in vivo. PYROXD1 oxidizes NAD(P)H to NAD(P)⁺ with tightly controlled kinetics, reducing one molecule of O₂ into H₂O₂ per turnover. PYROXD1 is expressed across a range of tissue types in humans, including skeletal muscle, and is conserved across many species.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of PYROXD1 antibody (31093-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-PYROXD1 transfected HeLa cells. This data was developed using the same antibody clone with 31093-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 31093-1-AP (PYROXD1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31093-1-PBS in a different storage buffer formulation.