

For Research Use Only

IGFBP7 Polyclonal antibody

Catalog Number: 31128-1-AP



Basic Information

Catalog Number:

31128-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0815

GenBank Accession Number:

NM_001553.2

GeneID (NCBI):

3490

UNIPROT ID:

Q16270-1

Full Name:

insulin-like growth factor binding protein 7

Calculated MW:

29 kDa

Observed MW:

29-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse kidney tissue, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31128-1-AP (IGFBP7 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.