

For Research Use Only

C18orf34 Polyclonal antibody

Catalog Number: 31148-1-AP



Basic Information

Catalog Number:

31148-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33847

GenBank Accession Number:

BC130535

GeneID (NCBI):

374864

UNIPROT ID:

Q5BJE1

Full Name:

chromosome 18 open reading frame

34

Observed MW:

95-102 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB: mouse spleen tissue, rat ovary tissue, rat spleen tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31148-1-AP (C18orf34 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.