

For Research Use Only

UBAP1L Polyclonal antibody

Catalog Number: 31189-1-AP



Basic Information

Catalog Number: 31189-1-AP	GenBank Accession Number: NM_001163692	Purification Method: Antigen affinity Purification
Concentration: 350 ug/ml	GeneID (NCBI): 390595	Recommended Dilutions: WB: 1:500-1:3000
Source: Rabbit	UNIPROT ID: F5GYI3	
Isotype: IgG	Full Name: similar to ubiquitin-associated protein 1 (predicted)	
Immunogen Catalog Number: AG34816	Calculated MW: 41 kDa	
	Observed MW: 41 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Caco-2 cells,
Species Specificity: human	

Background Information

UBAP1L is a fairly unknown protein. It is predicted to enable ubiquitin binding activity and is believed to be involved in the ubiquitin-dependent protein catabolic process through the multivesicular body sorting pathway. It is also predicted to be part of endosomal sorting complexes required for transport. UBAP1L contains a SOUBA domain, which begins at amino acid residue 269 and extends to the C terminus of the protein. This domain is similar to the one found in UBAP1, where it is known to interact with ubiquitin (PMID: 38293907).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Caco-2 cells were subjected to SDS PAGE followed by western blot with 31189-1-AP (UBAP1L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.