

For Research Use Only

C18orf32 Polyclonal antibody

Catalog Number: 31237-1-AP



Basic Information

Catalog Number:

31237-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35313

GenBank Accession Number:

NM_001199346.1

GeneID (NCBI):

497661

UNIPROT ID:

Q8TCD1

Full Name:

chromosome 18 open reading frame

32

Calculated MW:

9kDa, 76aa

Observed MW:

12 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : U2OS cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31237-1-AP (C18orf32 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.