

For Research Use Only

SCYL1 Polyclonal antibody, PBS Only

Catalog Number: 31302-1-PBS

Featured Product



Basic Information

Catalog Number:

31302-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35298

GenBank Accession Number:

NM_020680

GeneID (NCBI):

57410

UNIPROT ID:

Q96KG9

Full Name:

SCY1-like 1 (*S. cerevisiae*)

Calculated MW:

90 kDa

Observed MW:

110 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31302-1-AP (SCYL1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31302-1-PBS in a different storage buffer formulation.



WB result of SCYL1 antibody (31302-1-AP; 1:1200; incubated at room temperature for 1.5 hours) with sh-Control and sh-SCYL1 transfected HEK-293T cells. This data was developed using the same antibody clone with 31302-1-PBS in a different storage buffer formulation.