

For Research Use Only

VWA8 Polyclonal antibody, PBS Only

Catalog Number: 31348-1-PBS



Basic Information

Catalog Number: 31348-1-PBS	GenBank Accession Number: NM_015058.1	Purification Method: Antigen affinity purification
Size: 1 mg/ml	GeneID (NCBI): 23078	
Source: Rabbit	UNIPROT ID: A3KMH1	
Isotype: IgG	Full Name: KIAA0564	
Immunogen Catalog Number: AG35404	Calculated MW: 215kDa, 1905aa	
	Observed MW: 204-215 kDa, 117 kDa	

Applications

Tested Applications:
WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

VWA8, also known as KIAA0564, is a mitochondrial matrix-targeted protein having a putative ATPase activity (PMID: 34660594). It has 3 isoforms with the MW of 215, 204, and 117 kDa.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




HepG2 cells were subjected to SDS PAGE followed by western blot with 31348-1-AP (KIAA0564 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31348-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using VWA8 antibody (31348-1-AP) at dilution of 1:200 and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red). This data was developed using the same antibody clone with 31348-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 31348-1-AP (VWA8 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31348-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 31348-1-AP (VWA8 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31348-1-PBS in a different storage buffer formulation.