

For Research Use Only

NAGPA Polyclonal antibody

Catalog Number: 31365-1-AP



Basic Information

Catalog Number:

31365-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34860

GenBank Accession Number:

BC012194

GeneID (NCBI):

51172

UNIPROT ID:

Q9UK23

Full Name:

N-acetylglucosamine-1-phosphodiester alpha-N-acetylglucosaminidase

Calculated MW:

515aa,56 kDa; 309aa,34 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB: HeLa cells, HepG2 cells, L02 cells, U2OS cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31365-1-AP (NAGPA antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.