

For Research Use Only

# PRRC2C Polyclonal antibody

Catalog Number: 31383-1-AP



## Basic Information

Catalog Number:

31383-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35409

GenBank Accession Number:

BC006090

GeneID (NCBI):

23215

UNIPROT ID:

Q9Y520

Full Name:

BAT2 domain containing 1

Observed MW:

450 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:6000

IF/ICC: 1:200-1:800

## Applications

Tested Applications:

WB, IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

WB : A431 cells, HeLa cells, K-562 cells, PC-3 cells

IF/ICC : HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to Tris-Acetate gel system followed by western blot with 31383-1-AP (BAT2D1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using PRRC2C antibody (31383-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).