

For Research Use Only

# PRRC2C Polyclonal antibody, PBS Only

Catalog Number: 31383-1-PBS



## Basic Information

Catalog Number:

31383-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35409

GenBank Accession Number:

BC006090

GeneID (NCBI):

23215

UNIPROT ID:

Q9Y520

Full Name:

BAT2 domain containing 1

Observed MW:

450 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to Tris-Acetate gel system followed by western blot with 31383-1-AP (BAT2D1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31383-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using PRRC2C antibody (31383-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 31383-1-PBS in a different storage buffer formulation.