

For Research Use Only

INTS1 Polyclonal antibody

Catalog Number: 31428-1-AP



Basic Information

Catalog Number:

31428-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35428

GenBank Accession Number:

BC069262

GeneID (NCBI):

26173

UNIPROT ID:

Q8N201

Full Name:

integrator complex subunit 1

Observed MW:

250 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : HEK-293T cells, HeLa cells, Jurkat cells, NIH/3T3 cells

Background Information

INTS1 is a component of the Integrator (INT) complex, a complex involved in the small nuclear RNAs (snRNA) U1 and U2 transcription and in their 3'-box-dependent processing.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31428-1-AP (INTS1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.