

For Research Use Only

# Adiponectin Polyclonal antibody

Catalog Number: 31509-1-AP



## Basic Information

Catalog Number:

31509-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_144744.3

GeneID (NCBI):

246253

UNIPROT ID:

Q8K3R4

Full Name:

adiponectin, C1Q and collagen domain containing

Calculated MW:

26 kDa

Observed MW:

26 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

rat

Positive Controls:

WB : rat adipose tissue,

## Background Information

Adiponectin (AdipoQ), an adipocyte-derived hormone, is one of the most abundant adipokines in the blood circulation. Adiponectin modulates a number of metabolic processes, including improving insulin sensitivity and anti-inflammatory activity. The role of AdipoQ in reproduction is not yet fully understood, but the expression of AdipoQ in reproductive tissues has been observed in various animals and humans, including chicken testis, bovine ovary, and human placenta. Adiponectin exerts its effects by activating a range of different signaling molecules via binding to two transmembrane AdipoQ receptors, AdipoR1 and AdipoR2. AdipoR1 is expressed primarily in the skeletal muscle, whereas AdipoR2 is predominantly expressed in the liver. AdipoQ May play a role in cell growth, angiogenesis and tissue remodeling by binding and sequestering various growth factors.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Rat adipose tissues were subjected to SDS PAGE followed by western blot with 31509-1-AP (Adiponectin antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.