

For Research Use Only

LYRM5 Polyclonal antibody, PBS Only

Catalog Number: 31548-1-PBS



Basic Information

Catalog Number:

31548-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35889

GenBank Accession Number:

BC071769

GeneID (NCBI):

144363

UNIPROT ID:

Q6IPR1

Full Name:

LYR motif containing 5

Observed MW:

11 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

LYRM5, also known as ETFRF1, belongs to the Complex1_LYR-like Superfamily. LYRM5 is a mitochondrial protein and plays key roles in acetate metabolism, and in essential Fe-S cluster biogenesis. In vitro, LYRM5 was found to inhibit ETF by promoting the removal of flavin from the ETF holoenzyme, thus potentially regulating the rate of β -oxidation (PMID: 35402183, PMID: 36266680). The calculated molecular weight of LYRM5 is 11 kDa.

Storage

Storage:

Store at -80°C .

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-LYRM5 (IP:31548-1-AP, 4ug; Detection:31548-1-AP 1:500) with mouse heart tissue lysate 1840 ug. This data was developed using the same antibody clone with 31548-1-PBS in a different storage buffer formulation.