

For Research Use Only

C1orf123 Polyclonal antibody, PBS Only

Catalog Number: 31578-1-PBS

Featured Product



Basic Information

Catalog Number:

31578-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36156

GenBank Accession Number:

BC010908

GeneID (NCBI):

54987

UNIPROT ID:

Q9NWX4

Full Name:

chromosome 1 open reading frame

123

Observed MW:

25kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

C1orf123, also known as CZIB, is encoded by an open reading frame resulting from the splicing of 8 exons consisting of 160 amino acids. It has been identified as one of the human O-GlcNAc Transferase (OGT) interactors, suggesting that it is involved in post-translational modification (PMID: 24536041).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31578-1-AP (C1orf123 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31578-1-PBS in a different storage buffer formulation.



WB result of C1orf123 antibody (31578-1-AP; 1:12000; incubated at room temperature for 1.5 hours) with sh-Control and sh-C1orf123 transfected Jurkat cells. This data was developed using the same antibody clone with 31578-1-PBS in a different storage buffer formulation.