

For Research Use Only

# MXRA8 Polyclonal antibody

Catalog Number:31581-1-AP



## Basic Information

Catalog Number:

31581-1-AP

Concentration:

300 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33429

GenBank Accession Number:

BC017312

GeneID (NCBI):

54587

UNIPROT ID:

Q9BRK3

Full Name:

matrix-remodelling associated 8

Calculated MW:

49 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : COS-7 cells, mouse ovary tissue, rat ovary tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31581-1-AP (MXRA8 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.