

For Research Use Only

RGPD5 Polyclonal antibody

Catalog Number: 31636-1-AP



Basic Information

Catalog Number:

31636-1-AP

Concentration:

300 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36253

GenBank Accession Number:

NM_001037866

GeneID (NCBI):

84220

UNIPROT ID:

Q99666

Full Name:

RANBP2-like and GRIP domain containing 5

Calculated MW:

198kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

IHC: 1:50-1:500

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC: mouse liver tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 31636-1-AP (RGD5 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).