

For Research Use Only

# RGPD5 Polyclonal antibody, PBS Only

Catalog Number: 31636-1-PBS



## Basic Information

Catalog Number:

31636-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36253

GenBank Accession Number:

NM\_001037866

GeneID (NCBI):

84220

UNIPROT ID:

Q99666

Full Name:

RANBP2-like and GRIP domain containing 5

Calculated MW:

198kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 31636-1-AP (RGPD5 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31636-1-PBS in a different storage buffer formulation.