

For Research Use Only

CNTLN Polyclonal antibody, PBS Only

Catalog Number: 31698-1-PBS



Basic Information

Catalog Number:

31698-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36354

GenBank Accession Number:

GeneID (NCBI):

54875

UNIPROT ID:

Q9NXG0

Full Name:

centlein, centrosomal protein

Calculated MW:

162kDa

Observed MW:

162 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

CNTLN is also named C9orf101 and C9orf39. CNTLN is required for centrosome cohesion and recruitment of CEP68 to centrosomes. CNTLN is localized at the centrosome. The calculated molecular weight of CNTLN is 162 kDa.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 31698-1-AP (CNTLN antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31698-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 31698-1-AP (CNTLN antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31698-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CNTLN antibody (31698-1-AP) at dilution of 1:200 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002), CL594-Phalloidin (red). This data was developed using the same antibody clone with 31698-1-PBS in a different storage buffer formulation.