

For Research Use Only

LDAH Polyclonal antibody, PBS Only

Catalog Number: 31714-1-PBS



Basic Information

Catalog Number:

31714-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36372

GenBank Accession Number:

BC017473

GeneID (NCBI):

60526

UNIPROT ID:

Q9H6V9

Full Name:

chromosome 2 open reading frame 43

Observed MW:

32 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

LDAH, also known as C2orf43, is a newly identified LD protein. It has 4 isoforms and is abundantly expressed in tissues that predominantly store triacylglycerol (TAG), including liver and white and brown adipose tissue (PMID: 24357060).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31714-1-AP (C2orf43 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31714-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 31714-1-AP (LDAH antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31714-1-PBS in a different storage buffer formulation.