

For Research Use Only

# ADGB Polyclonal antibody, PBS Only

Catalog Number: 31756-1-PBS



## Basic Information

Catalog Number:

31756-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36246

GenBank Accession Number:

NM\_024694

GeneID (NCBI):

79747

UNIPROT ID:

Q8N7X0

Full Name:

chromosome 6 open reading frame  
103

Calculated MW:

189kDa

Observed MW:

235 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, IHC, IF-P, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

ADGB (Androglobin) is an evolutionarily conserved member of the globin protein family. It is specifically and highly expressed in spermatogenic cells (particularly sperm cells) of mammalian testes and in ciliated cells across various tissues. The ADGB protein possesses a unique structure, comprising a globin domain and a calmodulin-binding domain, suggesting its potential role in integrating gas (such as oxygen) sensing with calcium ion signaling functions. Current research indicates that ADGB plays a critical role in spermatogenesis within the male reproductive system, as well as in the formation and functional maintenance of cilia. Its functional deficiency may lead to abnormal sperm morphology and motility defects, thereby affecting male fertility. The target protein ADGB has a size of 235 kDa. This protein is unstable and prone to degradation, with degradation bands observed at 135 kDa, 70 kDa, and 55 kDa. (PMID: 36995441)

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 31756-1-AP (ADGB antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31756-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 31756-1-AP (C6orf103 antibody) at dilution of 1:100 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31756-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using 31756-1-AP (C6orf103 antibody) at dilution of 1:100 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31756-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using 31756-1-AP (C6orf103 antibody) at dilution of 1:100 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31756-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse testis tissue using ADGB antibody (31756-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31756-1-PBS in a different storage buffer formulation.