

For Research Use Only

KBTBD11 Polyclonal antibody

Catalog Number: 31778-1-AP



Basic Information

Catalog Number:

31778-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35389

GenBank Accession Number:

XM_011534772

GeneID (NCBI):

9920

UNIPROT ID:

O94819

Full Name:

kelch repeat and BTB (POZ) domain containing 11

Calculated MW:

65kDa

Observed MW:

68 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:2000-1:12000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Caco-2 cells, HEK-293T cells

Background Information

The putative tumor suppressor gene, KBTBD11, might play a role in tumorigenesis, and is associated with cellular apoptosis.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31778-1-AP (KBTBD11 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.