

For Research Use Only

ENGASE Polyclonal antibody

Catalog Number: 31784-1-AP



Basic Information

Catalog Number:	31784-1-AP	GenBank Accession Number:	NM_001042573	Purification Method:	Antigen affinity Purification
Size:	500 ug/ml	GeneID (NCBI):	64772	Recommended Dilutions:	IHC 1:50-1:500
Source:	Rabbit	UNIPROT ID:	Q8NFI3		IF/ICC 1:200-1:800
Isotype:	IgG	Full Name:	endo-beta-N-acetylglucosaminidase		
Immunogen Catalog Number:	AG34633	Calculated MW:	84 kDa		

Applications

Tested Applications:	IHC, IF/ICC, ELISA	Positive Controls:	
Species Specificity:	human, mouse	IHC :	mouse spleen tissue, mouse thymus tissue
		IF/ICC :	U-251 cells,

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Endo- β -N-acetylglucosaminidase (ENGASE) is an enzyme that belongs to glycoside hydrolase (GH) families 18 and 85, and is widely distributed in bacteria, fungi, and higher-order species. ENGASE hydrolyzes the β 1-4 linkage of the N, N-diacetylchitobiose (GlcNAc β 1-4GlcNAc) region of the N-linked glycans present on glycoproteins (PMID: 36435009). ENGASE is widely expressed, with higher expression in thymus and spleen (Refer to Uniprot).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse spleen tissue slide using 31784-1-AP (ENGASE antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Ethanol) fixed U-251 cells using ENGASE antibody (31784-1-AP) at dilution of 1:400 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002).