

For Research Use Only

CD134/OX40 Polyclonal antibody, PBS Only

Catalog Number: 31805-1-PBS



Basic Information

Catalog Number:

31805-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1261

GenBank Accession Number:

NM_003327

GeneID (NCBI):

7293

UNIPROT ID:

P43489

Full Name:

tumor necrosis factor receptor superfamily, member 4

Calculated MW:

29 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

CD134 (OX40) is a member of the TNFR-superfamily of receptors. Predominantly expressed on activated T cells, CD134 is activated by its cognate ligand CD134L (OX40L) and functions as a T cell co-stimulatory molecule. CD134-CD134L interactions have been proposed as a potential therapeutic target for treating autoimmunity (PMID: 26215166). CD134 can interact with TRAF2, TRAF3, and TRAF5.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MJ cells were subjected to SDS PAGE followed by western blot with 31805-1-AP (CD134/OX40 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31805-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 31805-1-AP (CD134/OX40 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 31805-1-PBS in a different storage buffer formulation.