

For Research Use Only

PABPC1L Polyclonal antibody, PBS Only

Catalog Number: 31907-1-PBS



Basic Information

Catalog Number:

31907-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36248

GenBank Accession Number:

NM_001124756

GeneID (NCBI):

80336

Full Name:

poly(A) binding protein, cytoplasmic

1-like

Calculated MW:

68 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

PAP1L, also known as PABPC1L, is a human gene that has been associated with various cellular processes, particularly in the context of reproduction and embryonic development. PABPC1L was markedly upregulated in RCC, and high PABPC1L expression correlated with unfavorable prognosis and resistance to ICB (PMID: 38382068).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31907-1-AP (PABPC1L antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31907-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 31907-1-AP (PABPC1L antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31907-1-PBS in a different storage buffer formulation.