

For Research Use Only

HDGFL1 Polyclonal antibody

Catalog Number: 31923-1-AP



Basic Information

Catalog Number:

31923-1-AP

Concentration:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36944

GenBank Accession Number:

NM_138574

GeneID (NCBI):

154150

UNIPROT ID:

Q5TGJ6

Full Name:

hepatoma derived growth factor-like

1

Calculated MW:

27 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human testis tissue,

Background Information

Hepatoma-derived growth factor-like 1 (HDGFL1, also known as HRP-1, PWWP1, and Hdgfrp1) is a member of the hepatoma-derived growth factor (HDGF) family related to tumorigenesis and the development of cancer (PMID: 18823380; 20509080). It is specifically expressed in the testis and activated in the nucleus (PMID: 22217008).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31923-1-AP (HDGFL1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.