

For Research Use Only

# PITHD1 Polyclonal antibody

Catalog Number:31973-1-AP



## Basic Information

Catalog Number:

31973-1-AP

Concentration:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37120

GenBank Accession Number:

BC017208

GeneID (NCBI):

57095

UNIPROT ID:

Q9GZP4

Full Name:

chromosome 1 open reading frame

128

Calculated MW:

24 kDa

Observed MW:

24 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, HL-60 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31973-1-AP (PITHD1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.