

For Research Use Only

# LANCL3 Polyclonal antibody

Catalog Number: 32087-1-AP



## Basic Information

Catalog Number:

32087-1-AP

Size:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34968

GenBank Accession Number:

BC093667

GeneID (NCBI):

347404

UNIPROT ID:

Q6ZV70

Full Name:

LanC lantibiotic synthetase component C-like 3 (bacterial)

Calculated MW:

46 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEL cells, TF-1 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32087-1-AP (LANCL3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.