

For Research Use Only

# SPATA5 Polyclonal antibody, PBS Only

Catalog Number:32104-1-PBS



## Basic Information

Catalog Number:

32104-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36838

GenBank Accession Number:

BC048217

GeneID (NCBI):

166378

UNIPROT ID:

Q8NB90

Full Name:

spermatogenesis associated 5

Calculated MW:

97, 86, 75 kDa

Observed MW:

100 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

## Background Information

AAA+ ATPases SPATA5 is an ATP-dependent chaperone part of the 55LCC heterohexameric ATPase complex which is chromatin-associated and promotes replisome proteostasis to maintain replication fork progression and genome stability (PubMed:38554706). SPATA5 plays an essential role in the cytoplasmic maturation steps of pre-60S ribosomal particles by promoting the release of shuttling protein RSL24D1/RLP24 from the pre-ribosomal particles (PubMed:35354024, PubMed:38554706).

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32104-1-AP (SPATA5 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32104-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 32104-1-AP (SPATA5 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 32104-1-PBS in a different storage buffer formulation.