

For Research Use Only

PRHOXNB Polyclonal antibody, PBS Only

Catalog Number:32118-1-PBS



Basic Information

Catalog Number:

32118-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37232

GenBank Accession Number:

NM_001105577

GeneID (NCBI):

646625

UNIPROT ID:

A6NGE7

Full Name:

parahox cluster neighbor

Calculated MW:

19kDa

Observed MW:

20 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

PRHOXNB may Catalyzes the stereoselective decarboxylation of 2-oxo-4-hydroxy-4-carboxy-5-ureidoimidazoline (OHCU) to (S)-allantoin. PRHOXNB is expressed at extremely low Levels.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32118-1-AP (PRHOXNB antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32118-1-PBS in a different storage buffer formulation.