

For Research Use Only

CD300c Polyclonal antibody, PBS Only

Catalog Number: 32141-1-PBS



Basic Information

Catalog Number:

32141-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC022279

GeneID (NCBI):

10871

UNIPROT ID:

Q08708

Full Name:

CD300c molecule

Calculated MW:

224 aa, 25 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Background Information

CD300c is a member of the CD300 family of receptors, which are type I transmembrane proteins expressed in both myeloid and lymphoid lineages and modulate immune responses by their stimulatory and inhibitory capabilities (PMID: 30485414). CD300c, also known as CLM-6, CMRF-35 and IgSF16, is composed of an IgV-like extracellular domain, a transmembrane domain and a short cytoplasmic tail. It is expressed on the surface of monocytes, macrophages, dendritic cells, neutrophils, some T and B lymphocytes (PMID: 1349532; 30498497). CD300c has been reported as an activating receptor expressed on monocytes and that it has a potential role in inflammatory responses (PMID: 23571507). It has also been identified as a T cell co-inhibitory molecule (PMID: 30498497).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 32141-1-AP (CD300c antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 32141-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 32141-1-AP (CD300c antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 32141-1-PBS in a different storage buffer formulation.