

For Research Use Only

# OSBPL6 Polyclonal antibody, PBS Only

Catalog Number: 32158-1-PBS



## Basic Information

Catalog Number:

32158-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37534

GenBank Accession Number:

BC052259

GeneID (NCBI):

114880

UNIPROT ID:

Q9BZF3

Full Name:

oxysterol binding protein-like 6

Calculated MW:

934 aa, 106 kDa

Observed MW:

106 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32158-1-AP (OSBPL6 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32158-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using 32158-1-AP (OSBPL6 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 32158-1-PBS in a different storage buffer formulation.