

For Research Use Only

# SORCS3 Polyclonal antibody, PBS Only

Catalog Number:32248-1-PBS



## Basic Information

Catalog Number:

32248-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37710

GenBank Accession Number:

NM\_014978.2

GeneID (NCBI):

22986

UNIPROT ID:

Q9UPU3

Full Name:

sortilin-related VPS10 domain  
containing receptor 3

Calculated MW:

136kDa, 1222aa

Observed MW:

160 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

Sortilin related VPS10 domain-containing receptor 3 (SorCS3) belongs to the vacuolar protein sorting 10 protein (VPS10p) receptor family. SorCS3 is a specific receptor for nerve growth factor (NGF) that regulates NGF signal transduction across membranes, and SorCS3 is mainly expressed in the hippocampus and cortex of the brain (PMID: 35393432).

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32248-1-AP (SORCS3 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32248-1-PBS in a different storage buffer formulation.