

For Research Use Only

# LINGO3 Polyclonal antibody

Catalog Number:32267-1-AP



## Basic Information

Catalog Number:

32267-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37289

GenBank Accession Number:

NM\_001101391

GeneID (NCBI):

645191

UNIPROT ID:

POC6S8

Full Name:

leucine rich repeat and Ig domain containing 3

Calculated MW:

64kDa

Observed MW:

67 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A549 cells, Calu-1 cells, U-251 cells, U-87 MG cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32267-1-AP (LINGO3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.