

For Research Use Only

LINGO3 Polyclonal antibody, PBS Only

Catalog Number:32267-1-PBS



Basic Information

Catalog Number:

32267-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37289

GenBank Accession Number:

NM_001101391

GeneID (NCBI):

645191

UNIPROT ID:

POC6S8

Full Name:

leucine rich repeat and Ig domain containing 3

Calculated MW:

64kDa

Observed MW:

67 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32267-1-AP (LINGO3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32267-1-PBS in a different storage buffer formulation.