

For Research Use Only

SEC23B Polyclonal antibody

Catalog Number: 32278-1-AP



Basic Information

Catalog Number:

32278-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35557

GenBank Accession Number:

BC005404

GeneID (NCBI):

10483

UNIPROT ID:

Q15437

Full Name:

Sec23 homolog B (S. cerevisiae)

Calculated MW:

86 kDa

Observed MW:

83-86 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells, K-562 cells, Raji cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32278-1-AP (SEC23B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.