

For Research Use Only

RIC1 Polyclonal antibody, PBS Only

Catalog Number: 32377-1-PBS



Basic Information

Catalog Number: 32377-1-PBS	GenBank Accession Number: BC136616	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 57589	
Source: Rabbit	UNIPROT ID: Q4ADV7	
Isotype: IgG	Full Name: KIAA1432	
Immunogen Catalog Number: AG36426	Observed MW: 150 kDa	

Applications

Tested Applications:
WB, IP, Indirect ELISA

Species Specificity:
human

Background Information

RIC1 (RAB6A GEF complex partner 1) is a protein that plays a critical role in intracellular membrane trafficking. It encodes a protein of 1056 amino acid (aa) residues, including a putative nuclear localisation sequence (PMID: 9099877). Ric1 is composed of two β -propellers and a C-terminal α -solenoid domain that remodels the Rab6 nucleotide-binding pocket to facilitate nucleotide exchange (PMID: 39632878). RIC1 proteins play important roles in vesicular trafficking, regulation of membrane dynamics, cell division and signalling, and are essential for maintaining normal cellular function.

Storage

Storage:
Store at -80°C .
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32377-1-AP (KIAA1432 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32377-1-PBS in a different storage buffer formulation.



IP result of anti-RIC1 (IP:32377-1-AP, 4ug; Detection:32377-1-AP 1:1000) with HeLa cells lysate 1880 ug. This data was developed using the same antibody clone with 32377-1-PBS in a different storage buffer formulation.